

circulationmaxx™ reviver

NEUROMUSCULAR STIMULATOR



Mains powered or battery powered
– Truly portable (at home or travelling)



High intensity output settings
for foot plates: 1 to 99



High intensity output settings
for TENS pads: 1 to 99



Programme timer: 1 to 60 mins



14 Electrical Muscle Stimulation
pattern cycles for foot plates



10 Electrical Muscle Stimulation
pattern cycles for TENS pads



Infrared, hand-held remote control



2 pairs of TENS electro pads

